

Inhalt für 6 ADOX CMS Filme • To develop 6 ADOX CMS films

ADOTECH II



B & W DEVELOPER for CMS 20 film

Inhalt: 100 ml Konzentrat für 1500 ml Arbeitslösung

Cont.: 100 ml conc. to make 1500 ml of working solution

Verdünnung: 1+14 = 33 ml + 467 ml Wasser für 2 Filme

Dilution: 1+14 = 33 ml + 467 ml of water for 2 films

Entwicklungszeit in Minutes:

Developing time in minutes:

ISO	20°C / 68°F
12/12'	8,5 Min.
20/14'	10 Min.



ADOX Fotowerke GmbH • 15526 Bad Saarow • Made in Germany • ADOX.DE

Chargen-Nr.: 130911